A message from Adam Morgan, M.D. (Laboratory Medical Director | SSM Health St. Mary's Hospital - Madison) ...

Ionized Calcium specimen type change and Blood Gas stability change Notification Date: August 5, 2021 Effective Date: August 24, 2021

Explanation: Specimen requirements and stability have changed due to new instrumentation.

Assay and Test ID	Current Specimen type and stability	New – Effective 8/24/2021 Specimen type and stability
	Serum SST     Stable 5 days refrigerated	Whole blood – Dark Green Lithium Heparin (no gel) or heparinized syringe
Ionized Calcium LAB01286		<ul> <li>Transport refrigerated if testing will be delayed &gt;30 minutes from collection</li> <li>Stable 1 day refrigerated</li> <li>Tube must be completely filled, uncentrifuged and unopened</li> </ul>
Arterial Blood Gas LAB01269	Whole blood – Dark Green Lithium Heparin (no gel) or heparinized syringe	Whole blood – Dark Green Lithium Heparin (no gel) or heparinized syringe
Venous Blood Gas		
LAB01889	Recommend: Transport specimen at room temperature and test within 30	Recommend: Transport specimen at room temperature and test within 30
Bicarbonate LAB01260	minutes of collection	minutes of collection
	Iced specimens accepted for 2 hours	Iced specimens accepted for <b>1 hour</b>
PH Venous LAB01582	after collection	after collection
PH CO2 Venous		
Mixed LAB01254		

If you have questions or concerns, please contact Adam Morgan M.D., Medical Director of Laboratory Services, at 608-258-6914, or Amy Little Soldier, Chemistry Technical Specialist, at 608-258-6910

Please communicate to colleagues and staff in your area accordingly.

Thank you, Adam Morgan, M.D. Laboratory Medical Director Pathology, SSM Health St. Mary's Hospital Madison 700 South Park Street Madison, WI 53715 Office: (608) 258-6914