

St. Joseph's/Candler Health System
Laboratory Point of Care Department
StatSensor Creatinine/GFR Meter Training Checklist

Trainee Name: _____

Date: _____

☐ Initial Training

☐ Retraining

Describe Meter

- ☐ **Highlight meter components and key features**
 - ☐ Color touch screen-illuminates
 - ☐ 30-sec analysis time
 - ☐ Corrects for Hct; up to 65
 - ☐ 1.2 µL sample volume; critical to clean and dry patient finger
 - ☐ Rechargeable batteries
- ☐ **Describe Docking Station**
 - ☐ 3 functions: upload, download, charge battery
 - ☐ Tote is not for meter storage
 - ☐ Describe lights (rt: battery charge, ctr: communication, lt: connected)
 - ☐ Show how to dock properly
- ☐ **Describe touch-screen and buttons**
 - ☐ On/off/ sleep
 - ☐ Show when buttons/touch screen are interchangeable
- ☐ **Describe strips**
 - ☐ Front fill, between the layers
 - ☐ Do NOT under-dose or re-dose; overdosing is not possible
 - ☐ Date vials when opened; 90-day expiration after opening
 - ☐ Refrigerate for long term storage
 - ☐ Keep vial closed
- ☐ **Demonstrate how to manipulate the meter**
 - ☐ Keep horizontal / keep strip well dry
 - ☐ Apply sample or QC until countdown starts
 - ☐ Position of barcode reader

Using the Meter

- ☐ **Introduce Quick Reference Guides**
- ☐ **Log on/Log off**
 - ☐ Speaker icon
 - ☐ Scanning methods; soft key, screen etc. / press and let go
 - ☐ Always check blue dialogue bar to see what meter is expecting next
 - ☐ 2-step logout
 - ☐ Automatic logout when docking / timed logout / logout required
- ☐ **QC**
 - ☐ QC Lockout
 - ☐ Values displayed or pass/fail
 - ☐ Expel first drop
 - ☐ Date QC vials when opened; 90-day expiration after opening
 - ☐ Explain reasons for failures/ how to follow-up (use comments & repeat)
- ☐ **Run a sample**
 - ☐ Check or enter strip lot
 - ☐ Scan patient barcode
 - ☐ Touch GFR to select; enter patient information appropriate for equation
 - ☐ Enter comments; *before* accepting results
- ☐ **Review Results**
- ☐ **Change Battery**
 - ☐ Update date and time
- ☐ **Cleaning**

Trainee Signature: _____

Date: _____

Trainer Signature: _____

Date: _____