Respiratory Specimen Collection Guide for Laboratory Testing

Covid- ED/IP/*STAT Outpatient collection

Nasopharyngeal (NP) or Mid-Turbinate Collection Storage: Room Temperature Place NP swab into dry sterile container. If using the blue conical tube, make sure to break the applicator shaft off at line. Send specimen to the Lab STAT Testing must be completed within 1 hour of

collection. Refrigerated: NA Frozen: NA

*Outpatient STAT collections at sites not performing Abbott ID Now will Route to ARUP with 24-96 hr. turn around time. If faster Turn around time is required, patient should be tested using POCT method or the patient sent to location where Abbott ID Now is available.

MRSA- Nasal

Liquid Amies- White top eSwab with included swab. Storage: Room Temperature transport to the lab as soon as possible (within 2 hours) Refrigerated 2-8°C for 48 hours



Point of Care Testing

Use only swab included in test kit. Follow test kit instructions for collection.

Covid- Routine Outpatient collection

ARUP Transport Media
With OP/Nasal swab (ARUP#56409)
With NP swab (ARUP #12884)
OR
NP, OP or Nasal Swab in UTM
Storage: Room Temperature or Refrigerated



2-30°C up to 2 days

(Transport to ARUP Frozen)



Influenza and RSV

Quidel/Copan Flexible mini-tip flocked Nasopharyngeal swab and UTM.

Storage: Room Temperature

RSV: 8 hours

Influenza A/B: 72 hours



FOR RSV TESTING

Patient 18-59 years will route to ARUP

Collect UTM

Storage/transport: Refrigerated 2-30°C

72 hours

Rapid Respiratory Panels

UTM and Nasopharyngeal swab

Storage: Refrigerated 2-8°C recommended following collection and while in transit if not tested immediately.

Updated Collection Requirement 1/7/2025 Rapid Strep A Throat- Swab Culture EZ Double with NO additive.

Collect using both swabs.

Room Temperature transport to the lab as soon as possible

Storage: Room Temp storage: 72 hours Refrigerated 2-8°C for 72 hours



Updated 12/16/2024

Contact Shannon Preece <u>preeces@slhs.org</u> with questions or concerns