

Microbiology Orderables Tip Sheet

Questions? Call Lab Customer Service at 615-875-5227 (5-LABS)

Abscess Culture

Potential Sources

- Drainage
- Pus (>1ml if possible)

Container

- Sterile screw top
- Syringe

Culture Bact Abscess

- LAB3087

Body Fluid Culture

Potential Sources

- Amniotic fluid
Container: Sterile screw top
- Pericardial fluid
Container: Bactec bottles
- Peritoneal fluid
Container: Bactec bottles
- Pleural fluid
Container: Bactec bottles

Culture Bact Body Fluid

- LAB269

Bronchial Culture

Potential Sources

- Bronchial lavage/wash
- TBNA

Container

- Sterile screw top

Culture Bact BAL

- LAB267

Deep Wound Culture (Subcutaneous Wounds)

Potential Sources

- Debridement
- Surgical site

Container

- Eswab

Culture Bact Wound Deep

- LAB503

Joint Culture

Potential Sources

- Synovial fluid or synovial tissue
 - Elbow
 - Hip
 - Knee
 - Shoulder
 - Wrist

Container

- Sterile screw top
- Syringe

Culture Bact Joint

- LAB6578

Respiratory/Sinus Culture:

C Resp

Potential Sources

- Sputum
- Tracheal aspirate
- Transtracheal scar

Container

- Sterile screw top

Culture Bact Respiratory

- LAB900

Respiratory/Sinus Culture

C Sinus

Potential Sources

- Maxillary sinus
- Sphenoid sinus

Container

- Sterile screw top
- Eswab

Culture Bact Sinus

- LAB3096

Staph. aureus Screen/ Vancomycin Resistant Enterococcus Screen Culture SA Screen

Potential Sources

- Nares
- Skin

Container

- Eswab

Culture Bact SA Screen

- LAB6590

Staph. aureus Screen/ Vancomycin Resistant Enterococcus Screen Culture VRE Screen

Potential Sources

- Rectal

Container

- Eswab

Culture Bact VRE Screen

- LAB6589

Superficial Wound Culture (Epidermis and Dermis Wounds)

Potential Sources

- Ear
- External eye
- Rash
- Skin
- Surgical site

Container

- Eswab

Culture Bact Wnd Spr

- LAB6576

Tissue/Surgical Culture

Potential Sources

- Biopsy
Container: Sterile screw top
- Bone
Container: Sterile screw top
- Surgical
Container: Eswab
- Tissue
Container: Sterile screw top

C Tissue and/or Surgical Culture

- LAB271

Yeast Culture

Potential Sources

- Genital
Container: Eswab
- Throat
Container: Eswab
- Urine
Container: Urine gray top

Culture from Throat, Urine, Genital

- LAB6579