

DRUG REFERENCE CHART

Key Facts about Critical Drugs

DRUG CLASS & DRUG	COMMON PRESCRIPTION BRAND AND NAME	DRUG METABOLITE	DETECTION TIMES	
Amphetamines/ Stimulants			Urine (Days)	Oral Fluid (Hours)
Amphetamine	Adderall®, Benzedrine®, Dexedrine®	Amphetamine	3-5	Up to 36
Methamphetamine	Desoxyn®, Vick's® Inhaler, Metabolite of Didrex®	Methamphetamine, Amphetamine (metabolite)	3-5	Up to 36
Methylphenidate	Ritalin®, Concerta®, Metadate, Methylin®, Focalin	Ritalinic Acid (metabolite)	Up to 2	Up to 36
Anticonvulsants			Urine (Days)	Oral Fluid (Hours)
Gabapentin	Gabaron®, Neurontin®	Gabapentin	1-2	Up to 36
Pregabalin	Lyrica®	Pregabalin	1-3	Up to 36
Barbiturates			Urine (Days)	Oral Fluid (Hours)
Amobarbital	Amytal, Tuinal	Amobarbital	4-6	Up to 36
Butalbital	Fiorinal®, Fioricet®	Butalbital	4-6	Up to 36
Pentobarbital	Nembutal®	Pentobarbital	4-6	Up to 36
Phenobarbital	Belladonna®, Luminal®	Phenobarbital	Up to 16	Up to 36
Secobarbital	Seconal®	Secobarbital	4-6	Up to 36
Benzodiazepines			Urine (Days)	Oral Fluid (Hours)
Alprazolam	Xanax®, Niravam, Xanor®, Tafil®, Alprolox®	Alpha-hydroxyalprazolam	2-4	Up to 36
Chlordiazepoxide	Librium®, Tropium®, Risolid®, Clipoxide®	Nordiazepam, Oxazepam, Demoxepam	2-4	Up to 36
Clonazepam	Klonopin®, Rivotril®	7-aminoclonazepam	2-4	Up to 36
Diazepam	Valium®, Apzepam®, Vival®, Apozeepam®	Oxazepam, Nordiazepam, Temazepam	2-4	Up to 36
Flunitrazepam	Rohypnol®, Fluscan®	7-aminoflunitrazepam	2-4	Up to 36
Flurazepam	Dalmadorm®, Dalmane®	Desalkylflurazepam, Hydroxyethylflurazepam	2-4	Up to 36
Lorazepam	Ativan®, Temesta®, Lorabenz®	Lorazepam	5-7	Up to 36
Medazepam	Nobrium®	Oxazepam	2-4	Up to 36
Midazolam	Versed®, Dormicum®, Hypnovel®, Flormidal®, Dormonid®	1-hydroxymidazolam	2-4	Up to 36
Oxazepam	Serax®, Seresta®, Serepax®, Sobril®, Vapen®, Opamox®, Oxascand®, Vaben®	Oxazepam	2-4	Up to 36
Temazepam	Restoril®, Normison®, Euhypnos®	Temazepam, Oxazepam	1-4	Up to 36
Triazolam	Halcion®, Novodorm®, Songar®	Alpha-hydroxytriazolam	2-4	Up to 36
Drugs of Abuse			Urine (Days)	Oral Fluid (Hours)
Bath Salts	N/A	Various		Up to 36
Cocaine	N/A	Benzoylcegonine (metabolite)	1-2	Up to 36
Heroin	N/A	6-monoacetylmorphine (metabolite), Morphine (metabolite)	1-3	Up to 36
Marijuana	Marinol®	Tetrahydrocannabinol (THC)	1-30	Up to 36
MDMA	N/A	MDMA, MDA (metabolite)	3-5	Up to 36

For more information, visit www.cordantsolutions.com, email info@cordanth.com, or call 1-855-895-8099

©2015 Cordant Health Solutions™. All Rights Reserved. Do not duplicate or distribute without written permission from Cordant Health Solutions.

DRUG CLASS & DRUG	COMMON PRESCRIPTION BRAND AND NAME	DRUG METABOLITE	DETECTION TIMES	
MDMA (Ecstasy)	N/A	MDMA, MOA (metabolite)	3-5	Up to 36
Nicotine	N/A	Cotinine		Up to 36
Phencyclidine (PCP)	N/A	Phencyclidine (PCP)	3-7	Up to 36
Spice	N/A	Various		Up to 36
Muscle Relaxants			Urine (Days)	Oral Fluid (Hours)
Carisoprodol	Soma®, Equanil®	Carisoprodol, Meprobamate (metabolite)	1	Up to 36
Cyclobenzaprine	Flexeril®, Lissiril	Cyclobenzaprine	3-8	Up to 36
Meprobamate	Soma, Miltown®	Meprobamate	2-3	Up to 36
Opioids			Urine (Days)	Oral Fluid (Hours)
Buprenorphine	Buprenex®, Suboxone®, Subutex®	Buprenorphine, Norbuprenorphine (metabolite)	2-4	Up to 36
Codeine	Tylenol® 3, Tylenol® 4, Tylenol® 5	Codeine, morphine (metabolite)	2-3	Up to 36
Fentanyl	Actiq®, Fentora™, Duragesic®	Fentinyll, Norfentanyl (metabolite)	1-2	Up to 36
Hydrocodone	Vicodin®, Lorcet®, Lortab®, Hycodan®, Norco®, Vicoprofen®, Anexsia®, Bekadid®, Hycodan®, Xodol®, Dicotid®	Hydrocodone, Hydromorphone (metabolite)	2-3	Up to 36
Hydromorphone	Dilaudid®, Exalgo®, Hymorphan®, Laudicon®, Novolaudin®, Palladone®	Hydromorphone	2-3	Up to 36
Meperidine	Demerol®, Algil®, Alodan®, Centrolgin®, Dispadol®, Dolantin®	Meperidine, Normeperidine (metabolite)	1-2	Up to 36
Methadone	Dolophine®, Methadose®	Methadone, EDOP (metabolite)	2-4	Up to 36
Naloxone	Narcan, Suboxone (combo)	Naloxone	2-4	Up to 36
Oxycodone	Tylox®, Percocet®, Percodan®, OxyContin®, OxyIR®, Roxicodone®	Oxycodone, oxymorphone	2-3	Up to 36
Oxymorphone	Opana®, Opana ER®, Numorphan®	Oxymorphone	2-3	Up to 36
Propoxyphene	Darvocet®, Darvon®	Propoxyphene, Norpropoxyphene (metabolite)	Up to 7	Up to 36
Tapentadol	Nucynta	Tapentadol		Up to 36
Tramadol	Ultram®, Ultracet®	Tramadol, O-desmethyltramadol (metabolite)		Up to 36
Other			Urine (Days)	Oral Fluid (Hours)
Alcohol	N/A	Ethyl Alcohol	Up to 24 Hours	Up to 36
Alcohol	N/A	Ethyl Glucuronide	80 Hours	Up to 36
Sleep Aids			Urine (Days)	Oral Fluid (Hours)
Zolpidem	Ambien®	Zolpidem	1	Up to 36
Zopiclone	Lunesta®	Zopiclone	1	Up to 36
Tricyclic Antidepressants			Urine (Days)	Oral Fluid (Hours)
Amitriptyline	Elavil®, Endep®, Tryptanol®, Etrafon®	Amitriptyline, Nortriptyline	4-10	Up to 36
Clomipramine	Anafrinil®	Norclomipramine	4-10	Up to 36
Doxepin	Adapin®, Prudoxin®, Sine-	Nordoxepin	2-8	Up to 36
Desipramine	Norpramin®, Pertofrane®	Desipramine	2-11	Up to 36
Imipramine	Tofranil®, Imiprin®, Deprinol®, Presamine	Imipramine, Desiripramine (metabolite)	2-5	Up to 36
Nortriptyline	Aventyl®, Pamelor®, Sensoval®, Nortrilen®, Allegron	Nortriptyline	4-19	Up to 36

NOTE: Detection times in Oral Fluid may extend out to two days depending on drug and extent of use.

Some benzodiazepines are detected into common metabolites such as Nordiazepam and Oxazepam. If a prescribed benzodiazepine is not listed, please consult our toxicologist.

References:

- Disposition of Toxic Drugs and Chemicals in Man. Baselt.
- Pain Management Testing Reference. White.

For more information, visit www.cordantsolutions.com, email info@cordanth.com, or call 1-855-895-8099

©2015 Cordant Health Solutions™. All Rights Reserved. Do not duplicate or distribute without written permission from Cordant Health Solutions.