

Specimen Collection Vials and Usages		
Vial	Product Name	Test Description(s)
	C&S Stool Vial Order #: 50388	<ul style="list-style-type: none"> <li>• Culture Stool Panel</li> <li>• Culture Stool Campylobacter</li> <li>• Culture Stool E. Coli</li> <li>• Culture Stool Salmonella &amp; Shigella</li> <li>• GI Panel PCR</li> <li>• PCR Enteric Bacterial Panel</li> </ul>
	C&S Urine Transport Kit Order #: 60313	<ul style="list-style-type: none"> <li>• Culture Urine Only</li> </ul>
	Clean Stool Vial Order #: 60378	<ul style="list-style-type: none"> <li>• Calprotectin (ARUP)</li> <li>• Clostridium difficile</li> <li>• Helicobacter pylori Antigen, Fecal (ARUP)</li> <li>• Lactoferrin</li> <li>• Norovirus</li> <li>• Rotavirus Antigen, Stool</li> <li>• Stool for WBCs</li> </ul>
	Ova & Parasite Kit Order #: 60381	<ul style="list-style-type: none"> <li>• Cryptosporidium Antigen, Feces</li> <li>• Giardia Antigen, Feces</li> <li>• Ova and Parasites Screen, Stool</li> </ul>
	H. Pylori Breath Kit	<ul style="list-style-type: none"> <li>• H. Pylori Breath Test (ARUP)</li> </ul>
	Aptima Urine Collection Kit Order #: 61739 <b>(Specimen must be filled within window on vial. If not, specimen is unacceptable)</b>	<ul style="list-style-type: none"> <li>• Chlamydia, Urine</li> <li>• Chlamydia/Neisseria, Urine</li> <li>• Neisseria, Urine</li> <li>• Trichomonas Urine (ARUP)</li> </ul>
	Urine Cup-Sterile Order #: 02068831	<ul style="list-style-type: none"> <li>• Chlamydia, Urine</li> <li>• Chlamydia/Neisseria, Urine</li> <li>• Neisseria, Urine</li> <li>• Trichomonas Urine (ARUP)</li> <li>• Routine or Complete Urinalysis</li> <li>• Urinalysis Micro</li> <li>• Culture Tissue</li> </ul>
	Thin Prep Pap Vial Order #: 60431	<ul style="list-style-type: none"> <li>• Liquid Based Pap</li> <li>• HPV HR and Genotyping 16/18/45</li> <li>• Chlamydia</li> <li>• Neisseria</li> <li>• Trichomonas</li> </ul>