

The Perfect Lab Specimen


Patient Name (Full First Name, Last name)
Ensure Correct Specimen Type
Requisition Complete
Full Birth Date DD/MM/YYYY
Each Specimen Labeled
Collector's Initials
Time and Date of Collection




Choose the correct tube for the test.



For All Liquid Specimens
 ✓ Tighten the screw tops before placing in specimen bags.
 ✓ Ensure that the bag is tightly closed.



Questions about Specimen Collection or Stability?
 Refer to the Online Test Directory.
<https://www.testmenu.com/concordhospitallab>



Label all specimens with the five requirements!

CONCORD HOSPITAL
Laboratory Services