
















Southwest General

Partnering with  University Hospitals

18697 Bagley Road, Middleburg Heights, OH 44130

LABORATORY SERVICES ORDER OF DRAW FOR MULTIPLE TUBE COLLECTIONS

CLOSURE COLOR	COLLECTION TUBE
	Blood Cultures - SPS
	Citrate Tube
 or 	Gel Separator Tube
	Serum Tube
	Serum Tube (RST)
 or 	Gel Separator Tube With Heparin
	Heparin Tube
 or 	EDTA Tube
	Separator Tube K ₂ EDTA with Gel
	Fluoride (glucose) Tube



Southwest General

Partnering with



University Hospitals

18697 Bagley Road, Middleburg Heights, OH 44130

LABORATORY SERVICES PROPER LABELING REQUIREMENTS

Primary Specimen Containers:

Specimen Containers are labeled with 2 patient-specific identifiers:

- Patient Name
- Date of Birth

Split samples must follow primary specimen requirements. Additional identifiers may include, but are not limited to:

- Patient Name
- Date of Birth
- Hospital #
- Social Security #
- Requisition #
- Accession # or unique random #

In limited situations, such as a trauma where a patient's identification is not known, a single identifier (Blood Bank ID # or MedTag #) may be used as long as it is traceable specifically to the trauma patient and uniquely identifies the specimen.

The Specimen identifier(s) must be indelible, legible and able to withstand all stages of processing and conditions of storage.