


Salem Health Laboratories

Order of Draw & Common Tubes

Order of Draw	Common Tubes	Additive	Number of Inversions
1		Sodium Polyanethole Sulfonate (SPS)	8-10
2		Sodium Citrate	3-4
3		Clot activator and gel (SST)	5
4		Clot activator (no gel)	5
5		Lithium Heparin and gel (PST)	8
6		Lithium Heparin or Sodium Heparin (no gel)	8
7		Spray Coated K ₂ EDTA	8
8		Spray Coated K ₂ EDTA	8
9		Sodium Fluoride/K Oxalate	8