

POUR OFF TUBES

Sterile Red Top Tube

Used for:

CSF
Any sterile body fluid



Blue Top Trace Metals Tube

Used for:

Trace element testing



Standard ARUP Tube

Used for:

Aliquoted specimen's that are being sent to ARUP



Protect from Light

Used for:

Any aliquoted specimen that needs to be protected from light.



WDL FROZEN ONLY TUBE
(Purple Top with Purple
writing)

Used for:

Any **FROZEN** WDL **BLOOD**
Specimen.



Standard Pour off Tube

Used for:

- WDL Refrigerated
- Mayo poured off specimens
- JC Virus to Quest
- Any in-house testing that needs to be poured off

Please try not to use the yellow caps on these tubes because the tubes can crack.

