


Microbiology Collection Device Guide

	<p>Cary Blair</p> <ul style="list-style-type: none"> • Stool Culture • Giardia Antigen 		<p>Total-Fix</p> <ul style="list-style-type: none"> • Ova & Parasite Screen
	<p>Sterile Cup</p> <ul style="list-style-type: none"> • Clostridium difficile Toxin A & B, PCR • Fecal Leukocytes • Sputum Culture • L. pneumophila, AG, EIA, Urine (LUA Antigen) 		<p>Hemocult Card</p> <ul style="list-style-type: none"> • Occult Blood, Stool <p>FOR ER/NICU Only the order is: POC Fecal Occult Blood Screen</p>
	<p>Aptima</p> <ul style="list-style-type: none"> • Chlamydia/GC Cervix, Endocervix, Rectal, Throat, Urethra, Vagina • Chlamydia/GC Urine Yellow Label Only 		<p>Sterile Tube with 0.5 cc Saline and Swab</p> <ul style="list-style-type: none"> • Wet Mount Trichomonas <p>FOR L&D Only: POC Wet Prep</p>
	<p>Red Top Swab</p> <ul style="list-style-type: none"> • MRSA PCR • SA NICU Surveillance • Strep A Screen w/RFX to Throat Culture • Aerobic Bacterial Culture (NOT FOR ANAEROBES) 		<p>Anaerobic Transport Device</p> <ul style="list-style-type: none"> • Anaerobic Culture, Aerobic, Fungal, and AFB Culture.
	<p>Universal Transport Media UTM (pink liquid/FLOQ Swab)</p> <p>Viral cultures to include HSV, Varicella, RSV, Influenza, Chlamydia, etc.</p> <p>***All of these are sent to a reference laboratory</p>		<p>Urine Culture Tube</p> <ul style="list-style-type: none"> • Urine Culture

Influenza, RSV, SARS COVID-19, RVP PCR

Viral Transport Media / Nasopharyngeal Swab



***Performed in house