





Microbiology Collection Device Guide

	Cary Blair - Stool Culture - Giardia Antigen		Total-Fix Ova & Parasite Screen
	Sterile Cup - Clostridium difficile Toxin A & B - Fecal Leukocytes - Sputum Culture L. pneumophila, AG, EIA, Urine (LUA Antigen)		Hemoccult Card Occult Blood, Stool FOR ER/NICU Only the order is: POC Fecal Occult Blood Screen
	Aptima Chlamydia/GC Cervix, Endocervix, Rectal, Throat, Urethra, Vagina Clamydia/GC Urine Yellow Label Only		Sterile Tube with 0.5 cc Saline and Swab - Influenza A&B Rapid Screen - Wet Mount Trichomonas FOR L&D Only: POC Wet Prep
	Red Top Swab - MRSA PCR - SA NICU Surveillance - Strep A Screen w/RFX to Throat Culture - Aerobic Bacterial Culture		Anaerobic Transport Device Anaerobic Culture If specimen collected in this device; it is also appropriate for fungal, AFB, and aerobic culture
	Universal Transport Media UTM (pink liquid/FLOQ Swab) Viral cultures to include HSV, Varicella, RSV, Influenza, Chlamydia, etc. All of these are sent to a reference laboratory		Urine Culture Tube Urine Culture