



Microbiology Collection Device Guide

	Cary Blair • Stool Culture • Giardia Antigen		Total-Fix • Ova & Parasite Screen
	Sterile Cup • Clostridium difficile Toxin A & B • Fecal Leukocytes • Sputum Culture • L. pneumophila, AG, EIA, Urine (LUA Antigen)		Hemocult Card • Occult Blood, Stool FOR ER/NICU Only the order is: POC Fecal Occult Blood Screen
	Aptima • Chlamydia/GC		Sterile Tube with 0.5 cc Saline and Swab • Influenza A&B Rapid Screen • Wet Mount Trichomonas FOR L&D Only: POC Wet Prep
	Red Top Swab • MRSA PCR • SA NICU Surveillance • Strep A Screen w/RFX to Throat Culture • Aerobic Bacterial Culture		Anaerobic Transport Device • Anaerobic Culture If specimen collected in this device; it is also appropriate for fungal, AFB, and aerobic culture
	Viral Transport (pink liquid) Viral cultures to include HSV, Varicella, RSV, Influenza, Chlamydia, etc. All of these are sent to a reference laboratory		Urine Culture Tube • Urine Culture