

## Using ARUP Gateway


### Browse A-Z

Click any letter in the search tool in the upper left corner of the window to display all tests beginning with that letter.



### Quick Search

Enter at least 3 characters of the Test Name/Number in the Search Box and select it from the Quick List.



### Full Search

Enter the Test Name or Number in the Search Box.



Click on the test to view the test information.

Results Found: 12	
Test Number	Test Name
81123	VITAMIN B2 - RIBOFLAVIN
92168	VITAMIN B3 (NIACIN)

### Viewing Test Information

Click on the tabs within the test to view additional information.



### Previously Viewed Tests

Click the History link to view as list of the last 50 tests viewed over the prior 50 days.



### Previously Viewed Tests

Test Number	Test Name	Last Viewed
565613	FOLATE, RBC	09/20/2013
0020056	Serum Viscosity	09/20/2013
HDL	HDL Cholesterol	09/20/2013
MIDH	IDH1/IDH2 Exon 4 Mutation	09/20/2013
265133	FLT3 MUTATION DETECTION BY PCR	09/20/2013

### Email or Print Test Information

You can email or print test information.

1. With the test information displaying, click either (for emailing) or (for printing).
2. Depending on your selection, enter or select the following:

For email, enter the relevant email message information in the **Send Test Link By Email** dialog and click **Send**.

For print, select the check boxes next to the information you want to print in the **Print Options** dialog and click **Print**.