



Penn Medicine

To: UPHS Physicians and Staff
From: Department of Pathology and Laboratory Medicine
Division of Precision and Computational Diagnostics, Molecular Pathology Laboratory
Vivianna Van Deerlin, M.D., Ph.D., Director
Christopher Watt, M.D., Ph.D., Associate Director
Salvatore Priore, M.D., Ph.D., Instructor

Date: November 18, 2019

Re: NEW ASSAY: *HTT* CAG Repeat Expansion Analysis for Huntington Disease

In November, CAG repeat expansion testing for the *huntingtin* gene (*HTT*) will be available in the Molecular Pathology Laboratory at the Hospital of the University of Pennsylvania for diagnostic, predictive and prodromal genetic testing for Huntington Disease, a neurodegenerative disease characterized by progressive motor, cognitive, and neuropsychiatric dysfunction. Huntington disease is caused by an expansion in the number of CAG repeats in exon 1 of *HTT*. The differential diagnosis for Huntington disease includes frontotemporal degeneration and amyotrophic lateral sclerosis. This PCR-based assay can detect and size alleles of up to 200 CAG repeats. Larger may be detected, but not accurately sized.

Result Reporting Categories

≤26 CAG repeats	Normal
27-35 CAG repeats	Intermediate allele
36-39 CAG repeats	Reduced penetrance allele
≥40 CAG repeats	Expansion associated with fully penetrant Huntington disease

Testing Information

Turnaround time: Up to 15 business days

Ordering: Huntington Disease (PxCode HTTCASE) in PennChart

Acceptable specimens:

- Peripheral blood collected in an EDTA (lavender) tube
- OraCollect-Dx oral swabs (Lawson # 281926)
- Frozen tissue (call the laboratory to discuss)

Informed Consent: A consent form, signed by the patient and the ordering physician/genetic counselor must be kept on file by the ordering physician/genetic counselor.

Contact Information: Call the Molecular Pathology Laboratory (215-662-6121) weekdays during regular business hours. For information on ordering and specimen requirements, kindly consult the [Lab Tests Services Guide](#).