



---

To: Outpatient Collection Facilities and Staff  
From: **Eline T. Luning Prak, MD, PhD**  
**Joyce Gonzalez, BS, M(ASCP)**

Date: 6/19/2023

Re: **New orderable for Urine Protein Electrophoresis with Reflex to Immunofixation**

---

In order to provide accurate results, an additional orderable has been added for Urine Protein Electrophoresis. Please order Urine Protein Electrophoresis, Random Urine for a random/spot urine. Please continue to order Urine Protein Electrophoresis Urine for a 24-hour collection.

---

UPEP, UIFE

### Collection Container

---



24 Hour Urine Container



Sterile, Leak-Proof Container

### Collect

---

Urine, 24 hr urine container with no preservative (NIL)

Also Acceptable: 100 ml random urine in a sterile, leak-proof container, but quantitation of protein peak will not be provided

For clinical questions related to this change, please contact the Immunology Resident at 215-980-9871. For operational questions related to this change, please contact Joyce Gonzalez at [Joyce.Gonzalez@pennterms.upenn.edu](mailto:Joyce.Gonzalez@pennterms.upenn.edu) or 215-662-6023.

[unsubscribe](#)