



**ESwab Regular Collection Kit (Lawson #151477)**

White cap w/Liquid Amies and a regular flocked swab

**Intended Uses:**

**All types of Anaerobic/Aerobic Cultures**

- |                            |  |
|----------------------------|--|
| <b>Fungal Cultures</b>     | <b>MRSA Nasal PCR</b>                              |
| <b>AFB Cultures</b>        | <b>SA/MRSA PCR (Nasal Preop)</b>                   |
| <b>Group B Strep PCR</b>   | <b>MRSA Wound PCR</b>                              |
| <b>Strep A PCR</b>         | <b>vanA VRE PCR</b>                                |
| <b>Mpox PCR</b>            | <b>Rotavirus *clean dry container is preferred</b> |
| <b>C. auris Screen PCR</b> |  |



**Nasopharyngeal (NP) UTM/VTM Media Collection Kit (Lawson #189185) 3 mL\***

Red cap with universal/viral transport media and flexible mini-tip flocked swab

- Intended Uses:**
- Respiratory Panel PCR w/ COVID 19**
  - COVID 19 Hospital PCR; COVID 19 BC PCR**
  - COVID 19-Flu A/B- RSV PCR, RSV/COVID PCR**
  - Bordetella pertussis/para PCR**
  - Legionella PCR**
  - Atypical pneumonia PCR (Mycopl/Chlamid)**
  - VZV PCR and HSV 1-2 PCR**
  - CMV QUAL PCR**

\*In case these tubes/swabs are not available, Nasopharyngeal (NP) swabs and UTM from other manufacturers are acceptable for all the above tests; swabs in Saline tubes are acceptable for all tests that include COVID.



**ESwab Flexible Mini-Tip Collection Kit (Lawson #156352)**

Blue cap w/Liquid Amies and a flexible mini-tip flocked swab

- Intended Uses:**
- Nasopharyngeal Culture**
  - Urethral Culture**
  - Bordetella pertussis/para PCR**
  - Legionella PCR**
  - Atypical pneumonia PCR (Mycopl/Chlam)**
  - VZV PCR and HSV 1-2 PCR**
  - CMV QUAL PCR**

While the blue swab is acceptable for some PCR tests the red swab is the swab of choice