








Test Name	Acceptable Collection	Information
<p>HPV High Risk</p>	 <p>Thin Prep Lawson# 158323</p>  <p>FLOQSwab for Self-Collection Lawson# 242100</p> <p>After collection, FLOQSwab® should be suspended in ThinPrep Pap Test PreservCyt® solution by trained personnel in clinical settings</p>	<p>HPV (Human Papilloma Virus) PCR test including HPV 16/18 Genotyping</p> <p><u>Co-testing</u> _____ 30-65 yrs of age</p> <p><u>ThinPrep PAP with HPV</u> Negative Cytology / HPV positive for HPV 16 or 18 can be referred directly to colposcopy</p> <p><u>Reflex for ASCUS</u></p> <p><u>ThinPrepPAP with ReflexHPV</u> ASCUS Cytology / HPV Positive (any High-Risk HPV type) can be referred directly to colposcopy</p> <p>Our HPV test has an Internal Control for sample adequacy.</p>
<p>HSV 1+2 PCR VZV PCR</p>	<p>Red cap e swab/UVT</p> 	<p>Lesion Swabs:</p> <p>HSV (Herpes Simplex Virus (1 and 2) and VZV DNA detection in Muco-cutaneous Lesions</p>
<p>Chlamydia/GC DNA (Outreach)</p>	<p>Roche Urine Kit Lawson# 176048</p> 	<p>Urine (Male and Female):</p> <p>Test for Chlamydia/GC, Trichomonas vaginalis, Mycoplasma genitalium</p>
<p>Trichomonas vaginalis PCR (Outreach)</p>	<p>Thin Prep Lawson# 158323</p> 	<p>ThinPrep:</p> <p>Test for Chlamydia/GC, Trichomonas vaginalis</p>
<p>Mycoplasma genitalium PCR (Outreach)</p>	<p>Roche Swab Kit Lawson#212521</p> 	<p>Genital, Vaginal, Endocervical Swab:</p> <p>Test for Chlamydia/GC, T. vaginalis, M. genitalium</p> <p>Anal, Throat Swab:</p> <p>Test for Chlamydia/GC</p>
<p>Vaginitis/ Vaginosis Panel PCR</p>	<p>Xpert Swab Kit Lawson# 169760</p> 	<p>Vaginal Swabs ONLY:</p> <ul style="list-style-type: none"> • Bacterial vaginosis • Candida group (<i>C. albicans</i>, <i>C. parapsilosis</i>, <i>C. tropicalis</i>, <i>C. dubliniensis</i>) • Candida glab-krus group (<i>C. glabrata</i>, <i>C. krusei</i>) • <i>Trichomonas vaginalis</i>