


Additional Herpes Simplex, Mycoplasma/Ureaplasma, Varicella Zoster, BV/CV Ordering Guide

 Sent to ARUP

Collection/Transport	Sample Type	Mycoplasma/ Ureaplasma by PCR LAB50106 ¹	Herpes Simplex by PCR LAB51182	Varicella Zoster by PCR LAB50097	BV/CV Wet Prep - 87210.08 ² Fungal Culture and KOH – 87102.06 ¹¹
ARUP UTM/VTM w/Swab	Genital Swab	X			
	Rectal Swab	X			
	Extragenital Vesicle		X	X	
	Upper Respiratory Swab	X			
ARUP UTM/VTM	Urine – 1 mL transferred to media	X			
Sterile Container or Transport Tube	Serum or EDTA Plasma		X	X	
	Ocular Fluid		X	X	
	CSF		X	X	
	BAL	X	X		
Rayon Culture Swab – Green Cap Swab – Red Cap Swab – Blue Cap					X

Notes:

1. Detects and speciates *Ureaplasma parvum*, *Ureaplasma urealyticum*, *Mycoplasma hominis*, and *Mycoplasma genitalium*; consider ordering for cases of non-gonococcal urethritis.
2. Screen for bacterial vaginosis, yeast, and *T. vaginalis* if present.
3. Culture and full identification of yeast, may assist in treatment decisions regarding fluconazole resistance.