





# Validated Specimen Matrix

Bellin

Collection Kit		Sample Type	Aptima Combo 2 (CT/NG)	Aptima Mycoplasma genitalium	Aptima Trichomonas vaginalis	Aptima HSV 1&2	Aptima BV	Aptima CV/TV	Aptima HPV	Aptima HPV Genotype
ThinPrep Pap Vial		Cervicovaginal	X		X				X	X <sup>1,2</sup>
Aptima Multitest Swab		Vaginal (Patient or Clinician Collected)	X		X					
		Penile Meatal (Patient Collected)								
		Throat	X							
		Rectal	X							
		Anogenital Lesion				X				
Aptima Unisex Swab		Endocervical	X <sup>3</sup>							
		Conjunctiva	X <sup>3</sup>							
		Male Urethral	X <sup>3</sup>							
Aptima Urine Specimen		Female Urine	X <sup>4</sup>							
		Male Urine	X <sup>4</sup>							

<sup>1</sup> Aptima HPV Genotype 16, 18/45 is automatically reflexed when pap negative and HPV positive, but can be added by request if both pap and HPV are positive

<sup>2</sup> Cytology department needs to be called when this test code is ordered

<sup>3</sup> Specimen sent to ARUP

<sup>4</sup> Patient should not have urinated at least 1hr prior to collection of urine

# Bellin Health Test Order Codes



Test Description	Test Code	Collection Device	Sample Type
Image-Guided Pap Test	88175.03	ThinPrep Pap Test	Cervicovaginal
Image-Guided Pap Test with reflex to Aptima HPV when ASCUS	88175.03	ThinPrep Pap Test	Cervicovaginal
Image-Guided Pap Test with reflex to Aptima HPV when ASCUS or Negative	88175.03	ThinPrep Pap Test	Cervicovaginal
Image-Guided Pap Test and Aptima HPV	87621.11	ThinPrep Pap Test	Cervicovaginal, anal/rectal
Aptima HPV Genotype 16, 18/45	LAB1107 <sup>1,2</sup>	ThinPrep Pap Test	Cervicovaginal, anal/rectal
Cytology Non-Gynecological	88112.04	ThinPrep Pap Test	Anal/rectal

Test Description	Test Code	Collection Device	Sample Type
Chlamydia trachomatis / Neisseria gonorrhoeae NAAT	LAB491	ThinPrep Pap Test	Cervicovaginal
	87491.11	Aptima Multitest Swab, Aptima Urine	Vaginal, Throat, Rectal, Urine
	LAB50232 <sup>3</sup>	Aptima Unisex Swab	Conjunctiva, Male Urethral, Endocervical
Trichomonas vaginalis NAAT	LAB1093	ThinPrep Pap Test	Cervicovaginal
	LAB1092	Aptima Multitest Swab	Vaginal
Chlamydia trachomatis / Neisseria gonorrhoeae / Trichomonas vaginalis NAAT	LAB1091	Aptima Multitest Swab	Vaginal
Herpes Simplex Virus 1 & 2 NAAT	LAB2097	Aptima Multitest Swab	Anogenital lesion

<sup>1</sup> Aptima HPV Genotype 16, 18/45 is automatically reflexed when pap negative and HPV positive, but can be added by request if both pap and HPV are positive

<sup>2</sup> Cytology department needs to be called when this test code is ordered

<sup>3</sup> Specimen sent to ARUP