

SWAB USAGE LIST FOR MICROBIOLOGY	
Swabs may be kept at room temperature 24-48 hours	
SWAB DESCRIPTION	TYPE OF CULTURE
E-Swab	Wound, Aerobic and Anaerobic
Double swaab with Sponge at bottom Red cap # 220109	Group B Strep culture GC Culture, female Nose Quick Strep - RGAS Rectal Throat Vaginal KOH Wound - Aerobic only Vaginal Culture
Mini-tip swab with Sponge at bottom Green cap #220133	GC Cultue, male Male genital (including GC) Stool Culture Ear Eye
Mini-tip (paper plastic wrapper)	Rapid Influenza Bordatella pertussis culture & PCR