

BD Liquid Amies Elution Swab (Eswab) Collection Transport System



SPECIFICATIONS:

STERILE: Yes

STORAGE/TRANSPORT TEMPERATURE

- 2 - 25°C

APPLICATION:

- For the collection of clinical specimens containing aerobes, anaerobes, and fastidious bacteria

COLLECTION

- Flocked swab in 1 mL Liquid Amies Elution

Aerobic/Anaerobic cultures acceptable for ESwab collection method:

- Wound/Abscess cultures
- Ear cultures
- Eye cultures
- Genital cultures
- Infection control samples:
Rectal/Axilla/nasal
- Fungal cultures – except body fluids including blood
- Tissue

Aerobic/Anaerobic cultures unacceptable for ESwab collection method:

- Any viral collection cultures
- Any Body Fluids:
 - CSF, Blood, PD, PT, Joint, etc.
- Streptococcus screen/Throat cultures
- Cystic Fibrosis cultures
- BAL cultures
- Stool cultures
- Trichomonas prep/antigen