

STOOL COLLECTION CHART

If the order is for

Stool Culture & Sensitivity
Bacterial Enteric Pathogens included:
Campylobacter
E. coli O157
Salmonella/Shigella
Yersinia



Para-Pak C & S (Culture & Sensitivity)
ORANGE CAP Vial (as shown above)
Room temperature – 24 hours
Refrigerated – 72 hours

If the order is for

Stool for Parasites
Amoeba
Cryptosporidium Antigen
Giardia Antigen
Giardia (not a duplicate of antigen)
Intestinal Parasites
Ova & Parasites
Worms



Para-Pak 10% Buffered Neutral Formalin
PINK CAP Vial (as shown above)
Room temperature – 24 hours
Refrigerated – 72 hours

If the order is for

***C.difficile* Toxin/Culture**
Fecal Lactoferrin
Rotavirus



Clean transport containers acceptable:
Clean transport container
Sterile urine cup
Refrigerate