

STOOL COLLECTION GUIDE

COLLECTION CONTAINER AND PEOPLESFT NUMBER	TEST NAME
Sterile Urine Cup PS# 15694 	• C. Difficile Toxin PCR (LAB0001372) Stability: 48 hours refrigerated • Helicobacter pylori antigen, feces (LAB0000325) Stability: 2 hours ambient, 72 hours refrigerated. Refrigerate within 2 hours of collection. • Calprotectin, Fecal (LAB0001972) Stability: 7 days refrigerated, 30 days frozen • Fecal Fat Qualitative (LAB0001273) Stability: 12 hours ambient, 2 weeks refrigerated
Stool Collection Container Obtain from Lab. 	• Fecal Fat Quantitative Collections 24 Hour- (LAB0003599) 48 Hour- (LAB0003600) 72 Hour-(LAB0003601) Stability: 12 hours ambient, 4 days refrigerated
IFOB PS# 72191 	• IFOB-Fecal Occult Blood Immunoassay (LAB0001314) Stability: 6 days ambient, 6 months refrigerated

COLLECTION CONTAINER AND PEOPLESOFT NUMBER	TEST NAME
Para-Pak Ultra EcoFix (Dark Green top) PS# 74499 	• Ova and Parasite (LAB000421) Stability: Ambient Indefinitely
ParaPak C & S (Orange Top) PS# 64670 	• Stool Culture, Gram Stain Optional (LAB00005533) Stability: 4 days ambient. Do Not Refrigerate. • Stool GI Panel (LAB0001619) Stability: 4 days ambient and refrigerated • C. Difficile Toxin PCR (LAB0001372) Stability: 48 hours refrigerated