










Lab Specimen Cheat Sheet		Anatomic pathology/ Surgical Pathology Lab				Core/Clinical Lab
		Intraoperative Evaluation (Frozen section/ gross margin evaluation/ lymphoma protocol)	Surgical Pathology Tissue Exam	Cytology <i>(Requires refrigeration if submitted afterhours)</i>	Leukemia/Lymphoma Panel (Flow Cytometry) <i>(Requires refrigeration if submitted afterhours)</i>	(Chemistry, microbiology, hematology, etc.)
Specimen	Tissue (ex. Biopsy, excision, resection, cores)	Fresh in specimen container 	Prefilled formalin jar or add formalin in pathology 	Not accepted	RPMI 	Fresh in sterile container 
	Fine needle aspiration, exfoliative brush (ex. TBNA, FNA, bronch brush, bile duct brush)	Not accepted	Not accepted	CytoLyt 	RPMI 	Fresh in sterile container 
	Body fluid (ex. BAL, washings, peritoneal/ pericardial/ pleural effusion)	Not accepted	Not accepted	Fresh in specimen container 	Fresh in specimen container 	Fresh in sterile container 