

Legacy Laboratory Services

Legacy LabAlert

February, 2024

An Important Update from Legacy Laboratory Services

BIOFIRE FILMARRAY Gastrointestinal Panel (GI PANEL) Norovirus

A recent alert from the manufacturer has reported a SLIGHT increase in false positive signals for the Norovirus target performed on the BIOFIRE GI panel.

A confirmation may be considered on samples detecting the Norovirus target, if medically necessary, by ordering the GeneXpert Norovirus PCR test.

If Norovirus is clinically suspected, it is recommended to order the GeneXpert Norovirus PCR assay to avoid delays in results.

Please see table 1 below for ordering guidance and specimen requirement.

Table 1:

EPIC Orderable	Cerner Orderable	Specimen Requirement
NOROV PCR	NOROV PCR	Random stool in a clean, unpreserved container

For additional information, please contact your account representative, client services or consult [the test table](#) on our website.

Legacy Laboratory Client Services: 503-413-1234, 877-270-5566, www.legacyhealth.org/labservices



LEGACY
LABORATORY
SERVICES

Save a copy to <https://legacyhealth.sharepoint.com/sites/LaboratoryServicesHub/LabAlerts/Forms/AllItems.aspx>