



# Supply Item No. 370167 – Sterile 4 oz Cup

Container	Label Placement Examples
<p>Preferred Patient Label placement:</p> <ul style="list-style-type: none"> <li>✓ Patient label applied over the manufacturer's barcode label</li> <li>✓ Applied vertically to container for ease of scanning</li> </ul>	
<p>NOT Preferred Patient Label placement:</p> <ul style="list-style-type: none"> <li>✗ Too close to manufacturer's barcode which can interfere with scanning the patient's barcode on the label</li> </ul>	
<p>Manufacturer's barcode label</p> <ul style="list-style-type: none"> <li>✗ Used by the manufacturer to track lots</li> <li>✗ Do not scan for patient testing</li> </ul>	