



Greiner Tubes

Chemistry Assays in **Mint** tube that can be performed on **Gold Top** during Tube Shortage

Albumin
ALK
ALT
Amylase
AST
Bili D
Bili T
Betahydroxy
BUN
Calcium
Chloride
Cholesterol
CO2 (Bicarb)
CPK
Creatinine
GGTP
Glucose
HDL
hsCRP
Iron
LDH
LDL m
Lipase
Magnesium
Phosphorus
Phenobarb
Potassium
Sodium
Total Protein
Transferrin
Triglycerides
Uric Acid

NO SUBSTITUTIONS!

MINT TUBE ONLY!

Troponin

NH3 - Ammonia

Please conserve Mint tops for these two tests!!

Substitution for Greiner Mint and Golds.



BD Gold