

Laboratory Technical Bulletin: Improved Intact PTH assay – Biotin Threshold Increase

Effective 10/6/25

Roche has updated their Intact PTH assay to increase its biotin threshold. The new assay display's no interference in samples with biotin concentrations up to 1200 ng/mL or 4920 nmol/L.

	Old	New
Assay measuring range	1.2 – 5000 pg/mL	6.0 – 5000 pg/mL
Assay Reference Range	14 – 65 pg/mL	17.9 – 58.6 pg/mL

Date submitted: September/30/2025

Submitted by: Ali Gabali MD PHD, Medical Director, McLaren Medical Laboratory
Victoria Coleman, System Manager Chemistry