

Microbiology Swab Collection Chart

BACTERIAL/FUNGAL COLLECTIONS



eSwab®

AFB culture, Anaerobe culture, Fungal culture, MRSA culture, Aerobic cultures except for cultures suspected or screening for GC



COPAN Dual Culture Swab System with Liquid Stuart's Medium


Eye/Throat screen for GC, Pre-Op Staph aureus MRSA Screen (SANPCR), Strep A (GASM), Routine Genital Culture for GC Screen

CHLAMYDIA AND GONORRHOEAE COLLECTIONS BY AMPLIFIED DETECTION




APTIMA Unisex Swab Specimen Collection Kits

Endocervical and Male Urethral Swab Specimens



APTIMA Multitest Swab Specimen Collection Kit

Vaginal Swab Specimens





VIRAL COLLECTIONS



Puritan Sterile Foam-Tipped Applicator Swab

Rapid Flu



Puritan HydraFlock Applicator Swab

Rapid RSV



Puritan brand Polyester-Tipped Applicator Swab

Viral Transport Media Collection Swab



Flocked Flexible Plastic Applicator Swab

Viral Transport Media NP Collection Swab



For more information regarding microbiology swab collections, please call (800) 934-6575 or visit NationwideChildrens.org/Lab.