Microbiology Swab Collection Chart

BACTERIAL/FUNGAL COLLECTIONS



eSwab® AFB culture, Anaerobe culture, Fungal culture, MRSA culture, Aerobic cultures except for cultures suspected or screening for GC



COPAN Dual Culture Swab System with Liquid Stuart's Medium Eye/Throat screen for GC, Pre-Op Staph aureus MRSA Screen (SANPCR), Rapid Strep A (GASM, RMGAS), Routine Genital Culture for GC Screen

CHLAMYDIA AND GONORRHOEAE COLLECTIONS BY AMPLIFIED DETECTION



APTIMA Unisex Swab Specimen Collection Kits Endocervical and Male Urethral Swab Specimens



APTIMA Multitest Swab Specimen Collection Kit Vaginal Swab Specimens



Viral Transport Media Collection Swab

Flocked Flexible Plastic Applicator Swab Viral Transport Media NP Collection Swab

For more information regarding microbiology swab collections, please call (800) 934-6575 or visit <u>NationwideChildrens.org/Lab</u>.

