





Urine Specimen Guidelines

Specimen Type	Testing	Stability		
		Urinalysis	Urine C&S	DAS
 Sterile Cup	Urinalysis Urine Culture & Sensitivities Drugs of Abuse Screen (DAS)	Room Temp 1 Hour Refrigerated 48 Hours Urine	Room Temp 2 Hours Refrigerated 24 Hours	Refrigerated 5 days
 Urine Transfer Tube (No Preservative)	Urinalysis Urine Culture & Sensitivities Drugs of Abuse Screen (DAS)	Room Temp 1 Hour Refrigerated 48 Hours Urine	Room Temp 2 Hours Refrigerated 24 Hours	Refrigerated 5 days
 Urinalysis Transfer Tube (with Preservative)	Urinalysis Only	Room Temp / Refrigerated 72 Hours	<i>Unacceptable Media</i>	<i>Unacceptable Media</i>
 Urine C&S Transfer Tube (With Preservative)	Urine Culture & Sensitivities Only	<i>Unacceptable Media</i>	Room Temp / Refrigerated 72 Hours	<i>Unacceptable Media</i>

For specimen collection questions, please call Avero Diagnostics Laboratory at 360.527.4580.

**focused on
answers.**

6221 Riverside Drive, Suite 119
 Irving, TX 75039
 877.232.9924

3560 Meridian Street, Suite 101
 Bellingham, WA 98225
 360.527.4580

averodx.com