


TUBE COLLECTION GUIDE- Venous Blood Collection

Additive		Additive	
	GOLD (SST Gel No Add)		WHITE (PPT -White Top EDTA)
Clot activator and gel for serum separation		K2EDTA and gel for plasma separation	
	RED (Red No Add)		GREEN (Gm Na H / •• Lt Gr Li)
Clot activator, Silicone coated		Sodium Heparin ••Lithium Heparin	
	LAVENDER (Lav EDTA)		PINK (Pink KEDTA)
Spray-coated K2EDTA		Spray-coated K2EDTA	
	LIGHT BLUE (BLUE NA C)		ROYAL BLUE (D Blu EDTA / •• D Blu None)
Sodium Citrate		Royal blue K2EDTA or Na2EDTA •• Royal blue (no additive or EDTA)	
	GRAY (Gry NaFl)		YELLOW ACD
Potassium Oxalate / Sodium Fluoride		ACD solution A / ACD solution B	
	VACUETTE 9NC COAG. (Blu/Wht Cit)		
Greiner VACUETTE® 9NC Coagulation Sodium Citrate 3.2%			

SPECIMEN LABELING PROCEDURE:

Please place the Epic label as shown in diagram below

- * All specimens must be labeled with an Epic generated label.
- * All specimens must be labeled immediately upon collection.



1. Patient's Name and Medical Record Number
2. Tube collection type
3. Collected date and time
4. Test(s) ordered