

URINE COLLECTION KIT

Order of Draw from BD Vacutainer® Urine Collection Cup

First Collection:

**BD Vacutainer® C&S Preservative Tube
(gray top)**

For Culture Only



Second Collection:

**BD Vacutainer® Conical UA Tube (yellow
& red speckled top)**

For Urinalysis Only

**Must be filled to 7.5 mL for proper urine/preservative ratio
If unsure whether enough urine is available to fill, use the clear top instead.**



Third Collection:

**BD Vacutainer® Hemoguard™ Tube PLH
(clear top, no additive)**

**For Urine Chemistry, Drug Screens, Pregnancy, UPEPs (2 tubes
needed), Low Volume Draws, Total Protein, Creatinine,
Microalbumin, etc.**

