
















	Current Stock Swab		New Stock Swab	Tests
BD Vaginal Self Collect		Aptima Multipurpose - Orange Vaginal Self Collect or Clinician Collect		PDC (chlamydia) PDG (gonorrhoeae) TRICD (Trichomonas)
BD Female Viper		Aptima Unisex - Purple Endocervical Male Urethral Rectal / Anal Clinician Collect Only		PDC (chlamydia) PDG (gonorrhoeae) TRICD (Trichomonas)
BD Male Viper		Aptima Unisex - Purple Endocervical Male Urethral Rectal / Anal Clinician Collect Only		PDC (chlamydia) PDG (gonorrhoeae) TRICD (Trichomonas)
		Aptima Urine		PDC (chlamydia) PDG (gonorrhoeae) TRICD (Trichomonas)
BD Male Urethral - mini tip		Aptima Urine		PDC (chlamydia) PDG (gonorrhoeae) TRICD (Trichomonas)
BD Blue Ames w/ charcoal - pertussis		eSwab - Mini-Tip – Blue nasopharynx urethral		AER (Bacterial Culture) ANA (Anaerobe Culture) BPERS (pertussis) GC (Gonorrhea Culture) FUN (Fungal Culture)
BD Tube Cultiurette Mini-Tip Green		eSwab - Mini-Tip – Blue nasopharynx urethral		AER (Bacterial Culture) ANA (Anaerobe Culture) BPERS (pertussis) GC (Gonorrhea Culture) FUN (Fungal Culture)

	Current Stock Swab		New Stock Swab	Tests
BD Red Dual Tip (Routine)				AER (Bacterial Culture) ANA (Anaerobe Culture) GC (Gonorrhea Culture) FUN (Fungal Culture) AFB (Mycobacterial) MRSA
BD White Single Tip - OB Routine		eSwab - Regular Collection - White throat, wound, etc.		
REMEL ACT 1 CULTURE TRANSPORT				

NO CHANGES IN SWABS LISTED BELOW				
	Affirm - Triple Screen			Vaginitis Vaginal Pathogens Triple Screen PVS
	BD Universal Transport Media			Influenza A and B FLRSV HSVRL VZVRL

For more information and complete collection details, visit our test menu:
www.testmenu.com/rochester