

















	Previous Stock Swab		New Stock Swab	Tests
BD Vaginal Self Collect		Aptima Multi-test - Orange Throat Rectal <b>Vaginal: Self Collect or Clinician Collect</b>		PDC (chlamydia) PDG (gonorrhoeae) TRICD (Trichomonas)
BD Female Viper		Aptima Unisex - Purple Endocervical Male Urethral Rectal / Anal <b>Clinician Collect Only</b>		PDC (chlamydia) PDG (gonorrhoeae) TRICD (Trichomonas)
BD Male Viper		Aptima Unisex - Purple Endocervical Male Urethral Rectal / Anal <b>Clinician Collect Only</b>		PDC (chlamydia) PDG (gonorrhoeae) TRICD (Trichomonas)
		Aptima Urine		PDC (chlamydia) PDG (gonorrhoeae) TRICD (Trichomonas)
BD Male Urethral - mini tip		eSwab - Mini-Tip – Blue nasopharynx urethral		AER (Bacterial Culture) ANA (Anaerobe Culture) BPERS (pertussis) GC (Gonorrhea Culture) FUN (Fungal Culture)
BD Blue Ames w/ charcoal - pertussis		eSwab - Mini-Tip – Blue nasopharynx urethral		AER (Bacterial Culture) ANA (Anaerobe Culture) BPERS (pertussis) GC (Gonorrhea Culture) FUN (Fungal Culture)
BD Tube Culturette Mini-Tip Green		eSwab - Mini-Tip – Blue nasopharynx urethral		AER (Bacterial Culture) ANA (Anaerobe Culture) BPERS (pertussis) GC (Gonorrhea Culture) FUN (Fungal Culture)
		Copan Flocked Swab		CMVQT (CMV Virus, neonatal saliva)



Previous Stock Swab	New Stock Swab	Tests
BD Red Dual Tip (Routine)		
BD White Single Tip - OB Routine		AER (Bacterial Culture) ANA (Anaerobe Culture) GC (Gonorrhea Culture) FUN (Fungal Culture) AFB (Mycobacterial) MRSA
REMEL ACT 1 CULTURE TRANSPORT		
 <p>eSwab - Regular Collection - White throat, wound, etc.</p>		

NO CHANGES IN SWABS LISTED BELOW		
Affirm - Triple Screen		Vaginitis Vaginal Pathogens Triple Screen (Candida, Gardnerella, Trich) PVS
BD Universal Transport Media		Influenza A and B FLRSV HSVRL VZVRL Gonorrhea Culture Chlamydia Culture

For more information and complete collection details, visit our test menu:  
[www.testmenu.com/rochester](http://www.testmenu.com/rochester)