

St. Joseph's/Candler Health System

Sure-Vue Urine Pregnancy Test Training Quiz

Employee Name (Please Print) _____ Signature _____

Operator ID _____ Department _____ Date _____

General Urine Pregnancy Testing Questions

1. Test kits are stored at what temperatures? _____ -- _____
2. Test cassettes are stable:
 - A. At Room Temperature
 - B. Until the printed expiration date on the box
 - C. If left in original pouch unopened
 - D. With passing QC after 30 days of storage
 - E. All of the above are true
3. Number the following steps (1-6) to correspond with the correct order in which they should be performed.

_____ Set a timer for three minutes

_____ Place the test device on a clean, level surface.

_____ Hold the dropper vertically and transfer urine sample

_____ Review and record results

_____ Allow test device & urine to equilibrate to room temperature

_____ Confirm patient identification on sample
4. How many drops of urine are used?
 - A. 2
 - B. 3
 - C. It doesn't matter as long as they the sample well is full
 - D. None of the above
5. Enter Determined Results:
 - A. If no bands appear in the test window, the results are: _____
 - B. If two bands appear in the test window, the results are: _____
 - C. If one band appears in the test window, the results are: _____

***To be completed by trainer/POC Representative**

POC Representative _____

Score ____ of ____ / ____ %