

St. Joseph's/Candler Point of Care
Clinitek Status Maintenance Log

Month/Year _____

Clinitek # _____

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
User Initials																															

Daily Maintenance:																															
Clean Exterior																															
Clean Test Table Insert																															

Weekly:																															
Disinfect Test Table																															
Inspect Calibration Bar																															

As Needed:																															
Clean Calibration Bar																															
Replace Printer Paper																															
Run Quality Controls																															

Document INU for days when the Clinitek Status is not in use.
Note: External Quality Controls must be run when opening a new box of cassettes and if the box is open 30 days.

Clinitek Status Maintenance Issue/Error Documentation

Date:Initials:Comments:

Monthly Review _____