

## Collection Tubes/Media



Gold Top/Red Marble Top with Gel Tube (SST) - for Serum Chemistry Assays



Plain Red Top Tube (No Gel/No SST) - for Serum Chemistry Assays



Light Blue Top Sodium Citrate Tube (2.7 mLs) - for Coagulation Assays



Fetal Fibronectin (FFN) Tube



24-Hour Urine Collection Container - for 24-Hour/Timed Urine Chemistry Assays



Liquid Amies (ESwab) - for Aerobic and Anaerobic Cultures



Lavender Top EDTA Tube (4 mLs) - for Routine Hematology



Pink Top EDTA Tube 6 mLs (Preferred) / Lavender Top EDTA Tube (10 mLs) - for Blood Bank Use



Light Blue Greiner Bio One Vacuette Sodium Citrate Tube (2.0 mLs) - for Platelet Function P2Y12-Plavix



Urine No Additive Tube - for Urinalysis



Culturette Swab - for Aerobic Cultures



Stool Culture Collection Vials



Light Green Top Lithium Heparin Tube (PST) - for Plasma Chemistry Assays



Dark Green Sodium Heparin Tube - for Plasma Chemistry Assays



Gray Top Sodium Fluoride Potassium Oxalate Tube - for Glucose Assay



Urine Preservative Tube - for Urine Culture



Anaerobic Culturette Swab - for Anaerobic Cultures



Stool O&P Collection Vials



Sterile Container



ThinPrep Vial



Viral Transport Medium



Cobas PCR Dual Swab Media - for GC/Chlam Probes (for Vaginal Specimens)



Aerobic (Blue Label) Anaerobic (Pink Label) Pediatric (Yellow Label)