Monday - Friday Urine (Catheterized): Day(s) / Time(s) 8am-5pm Performed: Turn Around Time\*\*: 24 hours **Specimen Collection Supplies:** • Specimen Collection Container If long delay is anticipated for delivery: • 50% Ethanol Optimum volume is 100ml. 50-100ml is adequate. Specimen Collection: Room temperature if specimen is delivered upon Specimen Handling: collection\* \*If overnight or weekend delay is anticipated, refrigerate \*If a longer delivery delay is unavoidable, add equal volume 50% ethanol Specimen Requirements: • Labeled with two (2) patient identifiers • Completed cytology specimen requisition • All specimens should be placed in a secondary "Biohazard Safety Bag". As Soon As Possible Deliver / Send to Teton **Transportation: Pathology** Catheterized urine is preferred for the lack of Comment / Special Instruction: contamination and greater content of transitional cells. If there are questions of low-grade transitional carcinoma, prior to insertion of catheter, it is helpful to collect voided baseline urine for comparison to the catheterized specimen. Excessive lubricant should be avoided as it may obscure the cells. Do not send first morning urine. **Rejection Criteria:** • Inadequate information / missing requisition • Unlabeled / Mislabeled Specimen • Specimen received in inappropriate fixative • First morning urine specimens

One (1) Week After Final Report is Issued

Cytopathology – Teton Pathology

Retention Time of Specimen:

**Laboratory Testing Performed:**