



# Ascension

## Ascension Michigan Laboratory Services

### Attention: ThinPrep Labeling Instructions

The new DYMO labels overlap on the bottom of the ThinPrep vials causing errors on the laboratory analyzer.

Beginning *IMMEDIATELY*:

- Labels are required to be placed HORIZONTALLY around the vial.
- Labels CANNOT cover the clear, blank portion of the vial.



*Horizontal Label*



*Clear, blank space unoccupied*

Thank you for your cooperation. Questions can be directed to the Cytology Department at 313-343-3011 or your Client Account Specialist.

19251 Mack Avenue  
Suite LL60  
Grosse Pointe Woods, MI 48236  
Phone: 800-863-5959  
Fax: 313-343-3854