BD Blood Collection Tube Conversion Crosswalk

Title: BD Blood Collection Tube Product conversion 2023

LOCATION: System **AUDIENCE**: All those collecting samples for blood testing

EFFECTIVE DATE: September 14, 2023 **APPROVED BY / DATE**: Shannon Preece 9/14/2023

PREPARED BY: Shannon Preece CONTACT: Shannon Preece

Program Manager for Specimen Management <u>preeces@slhs.org</u>

RATIONALE: St Lukes is converting by attrition to BD Vacutainer products for Blood Specimen Collection

Containers.

St Luke's

BD Tube Selection Old Greiner Tube Light Blue 3.0 mL Light Blue 2.7 mL (adult) Minimum Adult (black center) Lawson 122033 Fill volume required Light Blue 2.0 mL Pediatric Light Blue 1.8 mL (pediatric) Lawson 122034 (white center) Red plain 4.0 mL Red plain 4.0 mL Lawson 127450 Gold Top with gel 5.0 mL Gold Top 5.0 mL Lawson 202812 Green Lithium with Gel 3.5 mL Light Green Lithium with Gel 3 mL (yellow center) Lawson 202824 Dark Green Lithium Heparin3.0 mL Dark Green Lithium Heparin 4.0 mL Lawson 202821 (black center) Dark Green Sodium Heparin 4.0 mL Dark Green Sodium Heparin 4.0 mL (green center) Green Striped Label Lawson 202816 Lavender/EDTA 2.0 mL Lavender/EDTA 3 mL Lawson 118829 (white center) Pink/EDTA 6.0 mL Pink EDTA 6 mL Lawson 202818 Gray Na Fluoride Oxalate 2.0 mL Grey Top NA Fluoride/K Oxalate 4 mL (white center) Lawson 121822

- > Existing Greiner tubes should be used until supply is exhausted or expired.
- Refer to Epic Specimen Collection for number of tubes to draw.
- > Refer to St. Luke's Comprehensive Test Directory for draw order and instructions.