

Order of Draw

Greiner Tube Color	Tube Name	Mix by Inversion	BD Tube Color
	Blood Cultures	8 – 10 times *collect aerobic bottle first	
 *fill to black arrow	Light Blue (Sodium Citrate) Tubes	4 times	 *fill to minimum indicator line
	Red Or Gold (Serum) Tubes	5 – 10 times Clotting time: 30 minutes	
	Dark Green Or Light Green (Heparin) Tubes	8 – 10 times	
	Lavender Or Pink (EDTA) Tubes	8 – 10 times	
	Grey (Sodium Fluoride) Tube	8 – 10 times	