

## QUICK COMMUNICATION

TITLE: Chlamydia and GC collection changes			
LOCATION: SLHS Systemwide	AUDIENCE: All staff		
EFFECTIVE DATE: 9/18/2023	CONTACT: Obiageli Ozumba ozumbao@slhs.org Shannon Preece preeces@slhs.org		
PREPARED BY: Shannon Preece, Laboratory Program Manager	APPROVED BY / DATE: Obiageli Ozumba		

## Aptima Multitest Swab Specimen collection kit now validated/ FDA cleared for Throat and Rectal collection of CT/NG testing.

Thin Prep Media	Aptima Unisex Swab/Media	Aptima Urine Collection	Aptima Multitest Swab/Media
THE Prey Press And Press A	HOLOGIC* 29 rd. Aptimas* Grad Shortway Turneler Tubes	HOLOGIC 15 Aptima 20 mt.	HOLOGIC° 2.9 mL Aptima® Nagical Curils Transport Modes (STM)
<ul><li>Cervical</li><li>Endocervical</li></ul>	<ul><li>Endocervical</li><li>Male Urethral</li><li>Eye</li></ul>	• Urine	<ul><li>Vaginal</li><li>*Throat</li><li>*Rectal</li></ul>
	Lawson# 155790	Lawson #155789	Lawson# 199421 *New acceptable source for this swab type