

Title: Greiner Blood Collection Tube Product conversion October 2020

LOCATION: System	AUDIENCE: All those collecting samples for blood testing
-------------------------	---

EFFECTIVE DATE: October 13, 2020	APPROVED BY / DATE: Michele Harris, October 5, 2020
---	--

PREPARED BY: Michele Harris, Director Laboratory COE	CONTACT: Michele Harris, Director Laboratory COE MicheleH@slhs.org
---	---

RATIONALE: St Lukes is converting by attrition to Greiner Vacuette products for Blood Specimen Collection Containers. The Greiner descriptions and mSL label descriptions are below.

NEW GREINER Tube Selection

	Greiner Tube Description	mSL Label Tube Description
	<input type="checkbox"/> Light Blue 3.0 mL Adult (black center)	✓ Lt Blue 3.0mL –black cntr
	<input type="checkbox"/> Light Blue 2.0 mL Pediatric (white center)	✓ Lt Blue 2.0mL – white cntr
	<input type="checkbox"/> Red plain 4.0 mL	✓ Red 4.0mL
	<input type="checkbox"/> Gold Top with gel 5.0 mL	✓ Gold 5.0mL
	<input type="checkbox"/> Green Lithium with Gel 3.5 mL (yellow center)	✓ Grn Gel – yellow cntr
	<input type="checkbox"/> Dark Green Lithium Heparin 3.0 mL (black center)	✓ Grn Li Hep – black cntr
	<input type="checkbox"/> Dark Green <u>Sodium</u> Heparin 4.0 mL (green center)	✓ Grn Na Hep – green cntr
	<input type="checkbox"/> Lavender/EDTA 2.0 mL (white center)	✓ Lav/EDTA – white cntr
	<input type="checkbox"/> Pink/EDTA 6.0 mL	✓ Pink/EDTA 6.0mL
	<input type="checkbox"/> Gray Na Fluoride Oxalate 2.0 mL (white center)	✓ Gray 2.0mL – white cntr

*Tube containers will say ICE if specimen needs to be drawn on ice.

**Tube containers will say NO LITE if specimen needs to be protected from light.

*** SEE INSTRUCTIONS indicates that CTD needs to be consulted.

Title: Greiner Blood Collection Tube Product conversion August 2020

LOCATION: System

AUDIENCE: All those collecting samples for blood testing

EFFECTIVE DATE: August 11, 2020


APPROVED BY / DATE: Michele Harris, July 30, 2020

PREPARED BY: Michele Harris, Director Laboratory COE

CONTACT: Michele Harris, Director Laboratory COE
MicheleH@slhs.org

RATIONALE: As a result of a HealthTrust initiative, St Lukes is changing to the Greiner Vacuette products for Blood Specimen Collection Containers.

Old Tube Selection

	Light Blue 2.7 mL Light Blue 4.5 mL Light Blue 1.8 mL
	Red plain 6.0 mL Red plain 10.0 mL
	Gold Top 4.0 mL Tiger Top SST 8.5 mL
	Light Green Lithium with Gel 4.5 mL
	Dark Green Lithium Heparin 6.0 mL 2.0 mL
	Dark Green Sodium Heparin
	Lavender/EDTA 3 mL Lavender/EDTA 4 mL
	Pink EDTA 6 mL Pink EDTA 10 mL
	Grey Top – K Oxalate Grey Top – NA Fluoride

NEW GREINER Tube Selection

	Light Blue 3.0 mL Adult (black center)
	Light Blue 2.0 mL Pediatric (white center)
	Red plain 4.0 mL
	Gold Top with gel 5.0 mL
	Green Lithium with Gel 3.5 mL (yellow center)
	Dark Green Lithium Heparin 3.0 mL (black center)
	Dark Green <u>Sodium</u> Heparin 4.0 mL (green center)
	Lavender/EDTA 2.0 mL (white center)
	Pink/EDTA 6.0 mL
	Gray Na Fluoride Oxalate 2.0 mL (white center)

After the Greiner implementation the number of tubes available for use will be reduced;
 The number of tubes that need to be drawn will appear in specimen collection