

Vacurette Order of Draw

Tube Color	Tube Name	Mix by Inversion
	Blood Cultures	8 – 10 times *collect aerobic bottle first
	Light Blue (Sodium Citrate) Tubes	4 times *fill to black arrow
	Red Or Gold (Serum) Tubes	5 – 10 times
	Dark Green Or Light Green (Heparin) Tubes	8 – 10 times
	Lavender Or Pink (EDTA) Tubes	8 – 10 times
	Grey (Sodium Fluoride) Tube	8 – 10 times