



















UCDH MOLECULAR ACCEPTABLE SPECIMEN COLLECTION CONTAINERS

TEST	SPECIMEN COLLECTION CONTAINER	
<ul style="list-style-type: none"> CLL FISH Panel 	Sodium Heparin	
	EDTA	
<ul style="list-style-type: none"> HER-2 FISH (Breast & Gastric) 	FFPE Only	
<ul style="list-style-type: none"> Factor V Leiden Mutation Prothrombin (Factor II) DNA Mutation 	K2 EDTA	
<ul style="list-style-type: none"> Hereditary Hemochromatosis (HFE) Cystic Fibrosis DNA Mutations 	EDTA	
Viral Loads (Quantitative PCR) <ul style="list-style-type: none"> HIV-1 CMV HBV HCV EBV BK -Blood 	Plasma from EDTA	
<ul style="list-style-type: none"> BK Viral Load – Urine 	Cobas PCR Media – Urine Collection Kit	
<ul style="list-style-type: none"> Chlamydia Trachomatis Neisseria Gonorrhoeae (CT/NG) – URINE	Cobas PCR Media – Urine Collection Kit	

UCDH MOLECULAR ACCEPTABLE SPECIMEN COLLECTION CONTAINERS








TEST	SPECIMEN COLLECTION CONTAINER	
<ul style="list-style-type: none"> Chlamydia Trachomatis / Neisseria Gonorrhoeae (CT/NG) <p>FEMALE GENITAL – ENDOCERVICAL COLLECTION</p>	<p>Cobas PCR Media Dual Swab Kit ONLY (refer to kit instructions)</p>	
<ul style="list-style-type: none"> Chlamydia Trachomatis / Neisseria Gonorrhoeae (CT/NG) <p>FEMALE GENITAL – VAGINAL COLLECTION (Clinician and Self-Collections)</p>	<p>Cobas PCR Media Uni or Dual Swab Kit (Woven swab ONLY)</p>	
<ul style="list-style-type: none"> Chlamydia Trachomatis / Neisseria Gonorrhoeae (CT/NG) <p>ANORECTAL COLLECTION (Clinician and Self-Collections)</p>	<p>Cobas PCR Media Uni or Dual Swab Kit (Woven swab ONLY)</p>	
<ul style="list-style-type: none"> Chlamydia Trachomatis / Neisseria Gonorrhoeae (CT/NG) <p>THROAT COLLECTION (Clinician and Self-Collections)</p>	<p>Cobas PCR Media Uni or Dual Swab Kit (Woven swab ONLY)</p>	

UCDH MOLECULAR ACCEPTABLE SPECIMEN COLLECTION CONTAINERS

TEST	SPECIMEN COLLECTION CONTAINER	
<ul style="list-style-type: none"> Human Papillomavirus (HPV DNA) 	Hologic ThinPrep PAP TEST	<p><i>Important Notes:</i></p> <ol style="list-style-type: none"> Endocervical Brush and Spatula collections ONLY. Broom collection device is unacceptable for HPV DNA test on the Roche 4800 test system. ThinPrep PreservCyt Solution (NON-GYN) vial is unacceptable for HPV DNA test 
<ul style="list-style-type: none"> Respiratory Viral Panel 	BD UTM	
	M4	
	M4 RT	 
<ul style="list-style-type: none"> Bordetella Pertussis/Parapertussis DNA (continued on next page) 	BD UTM	

TO ENSURE ACCURACY OF MOST UP TO DATE INFORMATION, PLEASE DO NOT PRINT & DISTRIBUTE

UCDH MOLECULAR ACCEPTABLE SPECIMEN COLLECTION CONTAINERS

TEST	SPECIMEN COLLECTION CONTAINER	
<ul style="list-style-type: none"> Bordetella Pertussis/Parapertussis DNA (continued) 	M4	
	M4 RT	 
<ul style="list-style-type: none"> Covid-19 RNA Qual 	BD UTM	
	M4 RT	 
<ul style="list-style-type: none"> Covid-19 + Flu A/B 	BD UTM	
<p>For additional details regarding test requisitioning, specimen collection and transport requirements, results, and interpretation, search individual tests in our Lab Test Directory at https://www.testmenu.com/ucdavis</p> <p>Any further questions should be routed to Client Services Department: Phone: (916) 734-7373 Fax: (916) 737-1001</p>		