




## UCD DOPLM Tube Options for Transfusion Services

Test	Approved tube types for laboratory testing
<ul style="list-style-type: none"> <li>● Type and Screen</li> <li>● Crossmatch</li> <li>● Phenotype</li> <li>● A1 Lectin</li> <li>● Antibody Titer</li> </ul>	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>EDTA</p> </div> <div style="text-align: center;">  <p>RED TOP</p> </div> </div>
<ul style="list-style-type: none"> <li>● All Ortho Vision Testing (Automated Instrument Testing)</li> <li>● DAT</li> <li>● Fetal Maternal Hemorrhage (FMH)</li> <li>● Eluate</li> </ul>	<div style="text-align: center;">  <p>EDTA</p> </div>