





# UCD DOPLM HEMATOLOGY ALTERNATE SPECIMEN TUBES

TEST	Approved tube types for laboratory testing
<ul style="list-style-type: none"> <li>CBC W/ DIFF</li> <li>CBC W/O DIFF</li> <li>CBC W/ Man DIFF</li> <li>Hgb</li> <li>HCT</li> <li>Retic Studies</li> <li>ESR</li> <li>Rapid Malaria</li> <li>Blood Parasites</li> <li>Kleihauer Betke</li> <li>Hemoglobin Analysis</li> </ul>	 <p>EDTA      MAPP</p>
<ul style="list-style-type: none"> <li>PT/INR</li> <li>PTT</li> <li>DIC Screen</li> <li>Thrombin Time</li> <li>Factor Activities</li> <li>Inhibitor Screen</li> <li>von Willebrand testing</li> <li>Protein C and S testing</li> <li>Inhibitor Assay</li> <li>Anti-Xa testing</li> <li>Euglobulin Lysis Time</li> <li>PFA</li> <li>monomers</li> <li>TEG</li> <li>TEG Heparinase</li> <li>DOAC testing</li> <li>Fibrinogen</li> <li>D Dimer</li> </ul>	 <p>SODIUM CITRATE</p>
<ul style="list-style-type: none"> <li>TEG Platelet Mapping</li> </ul> <p>Both tubes required for platelet mapping tests to be performed</p>	 <p>SODIUM CITRATE      Lithium Heparin</p>
<ul style="list-style-type: none"> <li>Cryoglobulin</li> <li>HIT</li> </ul>	 <p>RED TOP</p>

## UCD DOPLM HEMATOLOGY ALTERNATE SPECIMEN TUBES

- FDP



Thrombin/Soybean  
Trypsin Inhibitor