





UCD DOPLM HEMATOLOGY ALTERNATE SPECIMEN TUBES

TEST	Approved tube types for laboratory testing
<ul style="list-style-type: none"> • CBC W/ DIFF • CBC W/O DIFF • CBC W/ Man DIFF • Hgb • HCT • Retic Studies • ESR • Rapid Malaria • Blood Parasites • Kleihauer Betke • Hemoglobin Analysis 	 <p style="text-align: center;">EDTA MAPP</p>
<ul style="list-style-type: none"> • PT/INR • PTT • DIC Screen • Thrombin Time • Factor Activities • Inhibitor Screen • von Willebrand testing • Protein C and S testing • Inhibitor Assay • Anti-Xa testing • Euglobulin Lysis Time • PFA • monomers 	<ul style="list-style-type: none"> • TEG • TEG Heparinase • DOAC testing • Fibrinogen • D Dimer  <p style="text-align: center;">SODIUM CITRATE</p>
<ul style="list-style-type: none"> • TEG Platelet Mapping <p>Both tubes required for platelet mapping tests to be performed</p>	 <p style="text-align: center;">SODIUM CITRATE Lithium Heparin</p>
<ul style="list-style-type: none"> • Cryoglobulin • HIT 	 <p style="text-align: center;">RED TOP</p>

UCD DOPLM HEMATOLOGY ALTERNATE SPECIMEN TUBES

- FDP



Thrombin/Soybean
Trypsin Inhibitor